




## Low Autofluorescence



### Conditions

[ Plate type ] InnoCell™ plate FP series ( flat bottom )  
 Non-treated ( N type )  
 [ Measurement ] Measured by Infinite® 200PRO M Plex

InnoCell™ plate has low autofluorescence at various wavelengths and even shows the values close to those of glass plates, which are known as an excellent material for fluorescence observation.

[ Abbreviation ] • PS : Polystyrene • FEP : Fluorinated Ethylene Propylene copolymer • PDMS : Poly ( dimethylsiloxane ) • COC : Cyclic Olefin Copolymer