





Low Autofluorescence

• Data obtained by Mitsui Chemicals


Excitation 360nm Blue fluorescence



Excitation 470nm Green fluorescence



Excitation 545nm Red fluorescence



Conditions

[Plate type] InnoCell™ T-plate FP series (flat bottom)
 Non-treated (N type)
 [Measurement] Plate reader

InnoCell™ T-plate has low autofluorescence at various wavelengths and clear fluorescent observation is also possible.

[Abbreviation] · PS : Polystyrene · PDMS : Poly (dimethylsiloxane) · FEP : Fluorinated ethylene-propylene